## SAVE THESE INSTRUCTIONS

- Never operate the unit if it has a damaged plug, if it is not working properly, if it has been dropped or damaged, or if it is immersed in water. Refer to an Authorized Service Center for repair service.
- DO NOT excessively bend or pull the cord. Unplug the cord by grasping the plug – NOT THE CORD.
- Keep the cord away from heated surfaces.
- Close supervision is necessary when this product is used by, on, or near children or disabled persons.
- Use this product only for its intended use as described in this manual. Do not use attachments unless recommended by the manufacturer.
- Never use while drowsy.

### WARNING

To reduce the risk of burns, electrocution, fire, or injury:

- DO NOT use the unit while bathing.
- ALWAYS unplug this product immediately after using.
- DO NOT place or store the unit where it can fall or be pulled into a tub, sink, or toilet.
- DO NOT reach for a unit that has fallen into water.
- Unplug immediately. Do not reach into the water under any circumstances. Have the unit checked by an authorized service center afterward.
- The transformer is designed for use in a 120 volt AC, 60 Hz outlet. Never force the plug into the outlet.
- DO NOT submerge the unit in water or any other liquid.

### ⚠ DANGER

To reduce the risk of electrocution:

## READ ALL INSTRUCTIONS BEFORE USING

## IMPORTANT SAFEGUARDS

- Turn Dental Water Jet to direct the jet of water toward tooth surfaces and between teeth to rinse away trapped plaque and particles. The pushing action of the water massages gums as it cleans the entire mouth.

**NOTE:** To prevent damage to delicate gum areas, do not point water jet directly at gums.

- Turn the unit off.

### Cleaning and Care

Thoroughly rinse reservoir with clean water after each use with dental rinse or mouthwash, and whenever the Dental Water Jet has not been used for an extended period of time. Allow full reservoir of lukewarm water to pass through power unit and cleaning tip before first use.

To purchase replacement tips, please contact our customer service. These are sold by pack of 2 tips under the reference WJT2RT.

## SAVE THESE INSTRUCTIONS

\*RBRC™ is a trademark of the Rechargeable Battery Recycling Corporation. The RBRC™ program provides a convenient alternative to placing used Ni-Cd batteries in the trash or the municipal waste stream, which may be illegal in your area. Please call 1-800-9-BATTERY™ for information on Ni-Cd battery recycling and disposal bas/re-strictions in your area. Conair's involvement in this program is part of our commitment to preserving our environment and conserving our natural resources."

Under various state and local laws, it is illegal to dispose of these batteries in our municipal waste stream. "The EPA certified RBRC® Battery Recycling Seal on the nickel-cadmium (NiCd) battery indicates Conair is voluntarily participating in an industry program to collect and recycle these batteries at the end of their useful life, when taken out of service in the United States or Canada. The RBRC® program provides a convenient alternative to placing used Ni-Cd batteries in the trash or the municipal waste stream, which may be illegal in your area. Please call 1-800-9-BATTERY™ for information on Ni-Cd battery recycling and disposal bas/re-strictions in your area. Conair's involvement in this program is part of our commitment to preserving our environment and conserving our natural resources."

### For Battery Disposal




- Contact your dentist prior to use if you have had oral or periodontal surgery in the previous twelve months.
- When you start any new oral hygiene program, your gums may initially bleed slightly. Contact your dentist if excessive bleeding occurs or if bleeding continues after two weeks of use.
- Be sure that only the water jet nozzle touches your teeth and gums. Damage to teeth and gums may occur if any other part of the unit touches them.
- This product is designed for oral cleaning only; DO NOT use the unit for any other purpose.
- As with any other new oral hygiene program, individuals with heart problems should check with their physician prior to using the Interplak® Dental Water Jet, and should always exercise particular care with any oral hygiene device.
- The Interplak® Dental Water Jet is a personal care unit and is intended for household use only.

### MEDICAL WARNING

For your safety:

### Instructions for Use

- Detach reservoir from power unit (see A).
- Open the jack cover on the power unit (see B).
- Plug the jack into the power unit and the transformer into a live outlet (see C) (make sure the power unit is turned off).
- Fully charge the unit for at least 16 hours (Note: If the outlet is controlled by a wall switch, be sure the switch is always on for continuous charge).
- Unplug the transformer from both the outlet and the power unit.
- Close the jack cover.
- Fit the water cleaning tip firmly into the power unit.
- Release and pull out the tube located at the base of the power unit (see D).
- Fill the reservoir with lukewarm water.
- Attach the reservoir to the bottom of the power unit, inserting the tube of the power unit into the reservoir, and aligning both arrows on the reservoir and the main unit (see E).
- Lock the reservoir into position, as indicated, by turning the reservoir until its arrow faces the locking position (see F).
- Turn the unit on by pressing the on/off switch.
- Hold mouth slightly open and elbow at mouth level to let water run back into the sink.

©2008 Conair Corporation, 150 Millford Road, East Windsor, NJ 08520  
For more information regarding this and other Interplak® products, visit [www.interplak.com](http://www.interplak.com)

Conair Corporation  
Service Department  
7475 N. Glen Harbor Blvd.  
Glendale, AZ 85307

East Windsor, NJ 08520

### Service Centers:

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

California residents need only provide proof of purchase and should call 1-800-334-4031 for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 12 months from the date of manufacture.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12 MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY.

Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

**Limited 1-Year Warranty**

Conair will repair or replace (at our option) your unit free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials. To obtain service under this warranty, return the defective product to the service center listed below that is nearest you, together with your purchase receipt and \$3.00 for postage and handling.

### Initial Preparation

Before using your Dental Water Jet for the first time :

- Carefully unpack the Dental Water Jet.
- Detach reservoir from power unit (see A).
- Open the jack cover on the power unit (see B).
- Plug the jack into the power unit (see C) and the transformer into a live outlet (make sure the power unit is turned off).
- Fully charge for at least 16 hours. (Note: If the outlet is controlled by a wall switch, be sure the switch is always on for continuous charge.)
- Unplug the transformer from both the outlet and the power unit.
- Close the jack cover.
- Fit the cleaning tip firmly into the power unit.
- Release and pull out the tube located at the base of the power unit (see D).
- Fill the reservoir with lukewarm water.
- Attach the reservoir to the bottom of the power unit, inserting the tube of the power unit into the reservoir, and aligning both arrows on the reservoir and the main unit (see E).
- Lock the reservoir into position, as indicated, by turning the reservoir until its arrow faces the locking position (see F).

- Turn the power unit on by pressing the on/off switch.
- Allow full reservoir of lukewarm water to pass through power unit and cleaning tip.
- Turn the unit off.

## Instruction Manual • Model WJ2B/WJ2CS

DENTAL WATER JET

# INTERPLAK®

by CONAIR

For your safety and continued enjoyment of this product, always read the instruction book carefully before using.

Questions or Comments? 1-800-334-4031  
Visit us at our website: [www.interplak.com](http://www.interplak.com)

071N7027  
IB-8131

## Welcome to a New Level of At-Home Dental Care

The Interplak® Water Jet enhances oral hygiene and significantly contributes to healthy gums when used as part of a responsible oral health management program. It is particularly effective at cleaning around orthodontic appliances, fixed and removable bridge work, implants and splints, and it provides valuable gum massage with every use.

Remember: Oral irrigation is a supplement to, not a substitute for, your oral hygiene routine. Ask your dental health care provider to advise you on how the Water Jet can most effectively meet your specific oral needs.

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